




SINGLE-USE MONOPOLAR FORCEPS

Continuous intraoperative nerve stimulation
integrated into surgical dissection.

- Enables real-time motor nerve mapping with triggered EMG
- Surgeon can evaluate nerve integrity without switching to a separate probe
- Common forceps patterns accommodate surgeon preferences and anatomical pathways
- Insulated instruments allows targeted, surgeon-directed stimulation
- Single-use design eliminates instrument reprocessing
- Compatible with all major intraoperative EMG systems
- Value-based pricing supports cost-effective use across procedures and departments



Scorpion combines ring forceps with continuous stimulation to support real-time motor nerve identification and functional assessment using evoked EMG. This approach reduces instrument exchanges by incorporating tEMG directly into surgical dissection. Multiple forceps patterns accommodate the surgeon's technique and instrument preference based on anatomy and surgical access.



NOW AVAILABLE IN SINGLE-USE

The Scorpion line offers four distinct ring forceps patterns for streamlined application of tEMG in complex procedures.

PROCEDURE-READY

Each Scorpion features an integrated lead wire and includes a stim return and ground needle for a complete, sterile single-use setup.

SELECT THE RIGHT PATTERN

Match the forceps pattern to the anatomy, surgical approach, and surgeon preference.

PATENTED COATING

Instruments feature a protective finish that provides electrical insulation and a durable, functional performance.

PROACTIVE® IONM, BUILT AROUND THE PROCEDURE

Together, Scorpion forceps, NMP procedure-based EMG electrodes, and the Nerveäna® Plus System create a PROACTIVE approach to IONM.

ORDERING INFORMATION

SINGLE-USE, SCORPION™ MONOPOLAR FORCEPS, BOX OF 5*

ITEM CODE	DESCRIPTION
I-DK-SHEM-M-5	Single-Use, Monopolar, Mosquito Forceps , Curved, Serrated, Fine Tip, 5" (12.7 cm)
I-DK-MCF512B-5	Single-Use, Monopolar, McCabe Forceps , Curved, Smooth, Fine Tip, 5.5" (14 cm)
I-DK-ATF714B-5	Single-Use, Monopolar, Adson Tonsil Forceps , Curved, Half-Serrated, Fine Tip, 7.25" (18.4 cm)
I-DK-MF512B-5	Single-Use, Monopolar, Mixer Forceps , Right-Angle, Serrated, Fine Tip, 5.5" (14 cm)

* Each single-use Scorpion includes a stim return and ground needle.

Scan for Scorpion Resources



tel: 866 815-6999
neurovisionmedical.com